OC SERIES OXYGEN CONCENTRATORS

EN

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
Käyttöohje

TABLE OF CONTENTS

Glossary.....3
 Special Notes4
 Before Using This Oxygen Concentrator4
 1. Introduction4
 2. Using Your OC Series Oxygen Concentrator8
 3. Cleaning, Maintenance and Servicing9
 4. Alarms and Troubleshooting Guide10
 5. Specifications.....11
 Appendix A: EMC Information12
 Appendix B: Circuit Diagram16
 Limited Warranty16

GLOSSARY

| | | | |
|---|------------------------------|---|---|
|  | Consult instructions for use |  | Follow Instruction for Use |
|  | Type BF Applied Part |  | Temperature limit |
|  | Power on |  | Power off |
|  | General warning Sign |  | Fragile, handle with care |
|  | Warning, electricity |  | AC Power |
|  | Stacking Limit By Number |  | Manufacturer |
|  | No open flames |  | Keep dry |
|  | No smoking |  | This Way Up |
| IP21 | Drip Proof Equipment |  | Class II symbol (Double insulated) |
|  | Humidity limit |  | Suitable for homecare use |
|  | Do not disassemble |  | Authorized representative in the European Community |
|  | Use no oil or grease | | |

SPECIAL NOTES

- Please read this manual carefully before using this product and save it for future reference.
- For assistance with this manual, contact your home care provider.
- Use the flow and duration settings of this product as prescribed by your physician.
- This equipment is not to be used as a life support device and supplies supplemental oxygen only.
- Children and patients with certain disabilities may be unable to understand or communicate product alarms, and therefore must be supervised when using this device.
- It is advised that patients maintain a backup oxygen supply in case of a machine malfunction or power shortage.
- If you experience an adverse reaction when using this device, contact a physician immediately.

BEFORE USING THIS OXYGEN CONCENTRATOR

- ⚠ **WARNING:** In high concentrations, oxygen vigorously accelerates combustion. Patients must make all efforts to reduce the risk of fire when undergoing oxygen therapy.
- ⚠ **WARNING:** This oxygen concentrator and all its components - including cannulas, connections and power cords - must be kept away from all sources of heat, open flame, sparks and static electricity.
- ⚠ **WARNING:** Do not smoke or allow others to smoke while you are undergoing oxygen therapy. Do not smoke or allow others to smoke in the same room as the oxygen concentrator and accessories.
- ⚠ **WARNING:** Oils and greases are prone to strong spontaneous combustion when exposed to oxygen under pressure. To avoid the risk of fire and personal injury, do not use oils and greases on or around the oxygen concentrator.
- ⚠ **WARNING:** This device is not suitable for use in the presence of a flammable anesthetic mixture with air, oxygen or nitrous oxide.

1. INTRODUCTION

Your physician has determined that supplemental oxygen is of benefit to you and has prescribed an oxygen concentrator set at a specific flow and duration setting to meet your needs. Do not change the flow or duration settings unless your physician tells you to do so.

Please read and understand this entire manual before using the device.

INTENDED USE

The OC Series Oxygen Concentrator is intended to provide supplemental oxygen to patients requiring oxygen therapy. The device is not intended to be life supporting or life sustaining.

ABOUT YOUR OC SERIES OXYGEN CONCENTRATOR

The OC Series Oxygen Concentrator delivers highly concentrated oxygen to patients requiring oxygen therapy. The device concentrates oxygen from the air for medical uses through a process called pressure swing absorption. Your home care provider will show you how to use your OC Series device and will be able to answer any questions you may have.

PARTS OF YOUR OC SERIES OXYGEN CONCENTRATOR

FIG 1: LCD PANEL



FIG 2: BUTTON BOARD

ACCESSORIES AND REPLACEMENT PARTS

Only use approved replacement parts.

The use of incompatible parts or accessories can result in reduced device performance. Contact your home care provider if you have any questions about the use of accessories or replacement parts.

WARNINGS AND CAUTIONS

A warning represents the possibility of harm to the operator or patient.

-  **WARNING:** For your concentrator to operate properly, air must be able to flow freely around the device. The ventilation ports are located at the rear base of the device and at the side air inlet filter. Keep the device at least 30 cm away from walls, furniture and especially curtains that may reduce airflow to the device. Do not place the concentrator in a small closed space (such as a closet).
-  **WARNING:** Do not remove the covers of this device. Servicing must be referred to an authorized and trained GCE home care provider.
-  **WARNING:** In the event of an equipment alarm or if you are experiencing any signs of discomfort consult your home care provider and/or your health care professional immediately.
-  **WARNING:** Oxygen generated by this concentrator is supplemental and should not be considered life supporting or life sustaining. In certain circumstances oxygen therapy can be hazardous; any user should seek medical advice prior to using this device.
-  **WARNING:** It is essential that any patient who is reliant on the oxygen supplied by this device has a backup source of oxygen in place.
-  **WARNING:** Do not use the oxygen concentrator if either the plug or power cord is damaged. Do not use extension cords or electrical adapters.
-  **WARNING:** Do not attempt to clean the concentrator while it is plugged into an electrical outlet.
-  **WARNING:** Do not connect the concentrator in parallel or in series with other oxygen concentrators or other oxygen therapy devices.
-  **WARNING:** Device operation above or outside of the voltage, flow, temperature, humidity and/or altitude values specified may decrease oxygen concentration levels.
-  **WARNING:** Your home care provider is responsible for performing appropriate maintenance at the intervals recommended by the device manufacturer.

CAUTIONS

A caution represents the possibility of damage to the equipment.

- Do not place liquids on or near the device.
- If liquid is spilled on the device, turn the power off and unplug from the electrical outlet before attempting to clean up the spill. Contact your home care provider if the device does not continue to work properly.

DANGER

To reduce the risk of burns, electrocution, fire or personal injury:

- Do not disassemble. There are no user-serviceable parts. Refer all servicing to an approved service professional.
- Avoid using the device while bathing. If your physician has prescribed continuous usage, the concentrator must be located in another at least 3 meters away from the bath.
- Do not come into contact with the concentrator while wet. Do not place or store the device where liquids can spill into the concentrator.
- Do not reach for the product if it has fallen into water. Unplug immediately and contact your product dealer.

RADIO FREQUENCY INTERFERENCE

This equipment has been tested and found to comply with EMC limited specified by IEC/EN 60601-1-2. These limits are designed to provide reasonable protection against electromagnetic interference in a typical medical installation.

TRANSPORT STORAGE AND OPERATING CONDITIONS

| | OPERATING | | TRANSPORT AND STORAGE | |
|----------------------|---|--------------------------------|---|-----------------------------------|
| TEMPERATURE |  | 10°C to 37°C (50°F to 98°F) |  | -30°C to 70°C (-22°F to 158°F) |
| RELATIVE HUMIDITY |  | 20% to 65%, no condensation |  | 15% to 95%, no condensation |
| ALTITUDE | 0 to 1828 meters | | - | |
| ATMOSPHERIC PRESSURE | 80 kPa - 101 kPa | | 50 kPa - 106 kPa | |
| ENVIRONMENT | Dry, well-ventilated, dust-free and pollution-free Away from electromagnetic interference | | Placed upright and vertical at all times | |

2. USING YOUR OXYGEN CONCENTRATOR

 **WARNING:** Do not use extension cords or electrical adapters.


1. Choose a location for your concentrator that allows it to draw in air freely. Make sure the device is at least 30 cm away from walls, furniture, curtains, or any other item that may restrict airflow to the device.
2. Install the concentrator's two stage filter:
 - a. Remove the filter access panel
 - b. Insert the inlet filter and turn clockwise to connect
 - c. Replace the filter access panel
3. Set up the humidifier bottle:
 - a. Press the humidifier connector button to remove the humidifier bottle
 - b. Remove the cap by turning the humidifier bottle clockwise
 - c. Fill the humidifier bottle with distilled or cooled boiled water between the min and max markings
 - d. Replace the humidifier cap and tighten until secure
 - e. Plug the humidifier bottle into the humidifier connector to reconnect it to the oxygen concentrator

NOTE: Replace the water in your humidifier bottle on a daily basis.

4. After reading and understanding the contents of this manual, plug the power cord into a grounded wall outlet and press the power switch to the On [I] position.
5. After turning on your OC-E Oxygen concentrator, allow at least 30 minutes for oxygen delivery to reach optimal concentration.
6. Make sure that the humidifier bottle is secure:
 - a. Use your finger to gently block the oxygen outlet on the humidifier bottle for 20 seconds
 - b. Once the float in the flow meter drops to the bottom of the gauge, remove your finger.
 - c. If the bottle makes a whistling sound - it means that the humidifier is properly secured to your device.
 - d. If you cannot hear this sound, remove the humidifier bottle, unscrew the cap, replace tightly and repeat the above test.
 - e. If no whistling sound is heard after taking these steps, contact your home care provider.
7. Adjust the oxygen flow so that the float mark is centered on the line marking the flow rate prescribed by your physician:

- a. Turn the flow control dial clockwise to increase the output flow rate.
 - b. Turn the flow control dial anti-clockwise to decrease the output flow rate.
8. Connect your oxygen tube to the oxygen outlet and put on your cannula as directed by your home care provider. If you are using an extension tube, connect this to the oxygen outlet and oxygen tube.
 9. When you are not using the oxygen concentrator, press the power switch to the Off [O] position to turn off the power.

OPERATION MODES AND LCD DISPLAY GUIDE

| MODE | GUIDANCE | LCD DISPLAY |
|------------------------|--|--|
| Start within 5 minutes | In the first five minutes of operation, the oxygen purity indicator, hour, minute and separator flicker at the same time. |  <p>Model OC-S has no oxygen indicator Model OC-E</p> |
| Continuous operation | The LCD indicates elapsed time and current run time. |  |
| Timed setting | <p>Set the length of your treatment by pressing the up and down buttons on the button board to increase or decrease the duration in 10-minute increments.</p> <p>20 seconds after setting your treatment duration, the timer will automatically begin to count down to the end of your treatment. The oxygen flow will stop once the timer runs out.</p> |  |

3. CLEANING, MAINTENANCE AND SERVICING

Always press the power switch to the Off [O] position and unplug your OC Series Oxygen Concentrator before cleaning.

The exterior of your oxygen concentrator should be cleaned weekly using a damp cloth with a mild household cleaner. Always wipe the exterior dry after cleaning.

Your cannula and humidifier must be cleaned and replaced as instructed by the manufacturer or your home oxygen provider.

You must clean the air filter in the OC Series Oxygen Concentrator at least once a week.

1. Open the access panel to the air filter to the side of the oxygen concentrator
2. Remove the air filters from the device
3. Rinse and allow to dry naturally
4. Once dry, reinstall the air filters

Maintenance and Servicing

The device has no user-serviceable parts, and must only be serviced by an authorized provider.

4. ALARM AND TROUBLESHOOTING GUIDE

The guide below lists alarms and common problems, as well as why they happened and what you can do to resolve them. If you are unable to resolve a problem or if you experience an issue that is not listed below, please contact your home care provider.

ALARMS

| ALARM | WHY IT HAPPENED | WHAT TO DO |
|--|--|--|
| Malfunction code H:07 is shown on the LCD display. | Device pressure is higher than specification. | Contact your home care provider. |
| Malfunction code H:06 is shown on the LCD display. | Device pressure is lower than specification. | Contact your home care provider. |
| Malfunction code H:05 is shown on the LCD display. | The device has overheated. | Contact your home care provider. |
| Malfunction code E:01 is shown on the LCD display | Output flow is too low. | Contact your home care provider. |
| Malfunction code E:02 is shown on the LCD display. | Thermistor failure. | Contact your home care provider. |
| Malfunction code E:03 is shown on the LCD display. | OSD sensor circuit failure. | Contact your home care provider. |
| Power failure alarm. | The device has detected a power outage or electrical supply failure. | Connect to a backup oxygen source and contact your home care provider. |

TROUBLESHOOTING GUIDE

The guide below lists common problems, as well as why they happened and what you can do to resolve them. If you are unable to resolve a problem, or if you experience an issue that is not listed below, please contact your home care provider.

| PROBLEM | WHY IT HAPPENED | WHAT TO DO |
|--|--|---|
| The device is not working when it is turned on, but the LCD display is showing. | Internal part failure. | Contact your home care provider. |
| The device is not working when it is turned on and the LCD display is not showing or appears intermittently. | The power cord is damaged. | Check whether the power cord is damaged. |
| | The power cord plug is not properly inserted into the electrical outlet. | Make sure the device is properly plugged in to the electrical outlet. |
| | The oxygen concentrator is not receiving power from the electrical outlet. | Check your household outlet fuse or circuit. |
| | Internal part failure. | Contact your home care provider. |

| | | |
|--|--|---|
| Oxygen is not flowing or the oxygen flow is limited. | The oxygen tubing is kinked or blocked, stopping the delivery of oxygen. | Unkink the oxygen tubing or replace tubing if necessary. |
| | The humidifier bottle is not properly connected to the device. | Reinstall the humidifier bottle. |
| The float in the flow meter does not move up or down when adjusting the flow control dial. | The flow control dial is not open. | Turn the flow control dial slowly and carefully. |
| | The flow control dial is faulty. | Contact your home care provider. |
| There is water in the cannula. | There has been a change in temperature; OR | Dry the inside of the humidifier cover. |
| | | Do not use hot water in the humidifier bottle. |
| | The device is too close to a wall, curtains or furniture. | Do not overfill the humidifier bottle. |
| | | Keep the oxygen concentrator and cannula in the same room at the same temperature. If using an extension hose with your oxygen concentrator, speak to your home care provider about fitting a water trap to collect excess moisture. |

5. SPECIFICATIONS

Technical Specifications

| Model | OC-E80 | OC-E100 |
|--------------------------|---|--|
| Power consumption (W) | 650 750W 6A (North America) | 650 750W 6A (North America) |
| Rating current (A) | 2.25 | 2.25 |
| Electrical requirements | 230V (+10%, - 15%) 50Hz (±1Hz) 230V (+10%, - 15%) 60Hz (±1Hz) 120V (+10%, - 15%) 60Hz (±1Hz) | 230V (+10%, - 15%) 50Hz (±1Hz) 230V (+10%, - 15%) 60Hz (±1Hz) 120V (+10%, - 15%) 60Hz (±1Hz) |
| Flow rate (l/min) | 0,5 to 8 | 1 to 10 |
| Purity | 87 – 96% | 87 – 96% |
| Outlet pressure (mpa) | 0.05±0.005 | 0.05±0.005 |
| Dimensions | 381 mm x 347 mm x 689 mm | 381 mm x 347 mm x 689 mm |
| Net weight | 27 kg | 27 kg |
| Sound level | 52dB(A) | 52dB(A) |
| Equipment class and type | Class IIa | Class IIa |
| EMC compliance | IEC 60601-1-2 | IEC 60601-1-2 |
| Features | Standard: LCD display; oxygen concentration indicator; temperature alarm; high/low pressure alarm; low oxygen concentration alarm; power failure alarm; timed setting; intelligent fault diagnosis; maintenance reminder. Additional: Positive pressure outlet | |

Disposal

Dispose of this device in accordance with local regulations.

APPENDIX A: EMC INFORMATION

If you have any questions regarding the guidance and declarations listed below, please contact your home care provider.

Guidance and Manufacturer's Declaration - Electromagnetic Emissions: This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure that it is used in such an environment.

| EMISSIONS TEST | COMPLIANCE | ELECTROMAGNETIC ENVIRONMENT-GUIDANCE |
|--|------------|---|
| RF emissions CISPR 11 | Group 1 | This device uses RF energy for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. |
| RF emissions CISPR 11 | Class B | The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network. |
| Harmonic emissions IEC 61000-3-2 | Class A | |
| Voltage fluctuations/Flicker emissions IEC 61000-3-3 | Complies | |

Guidance and Manufacturer’s Declaration - Electromagnetic Immunity: This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure that it is used in such an environment.

| IMMUNITY TEST | IEC 60601 TEST LEVEL | COMPLIANCE LEVEL | ELECTROMAGNETIC ENVIRONMENT-GUIDANCE |
|--|--|---|---|
| Electrostatic Discharge (ESD) IEC 61000-4-2 | ±8kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air | ±8kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%. |
| Electrical Fast Transient/Burst IEC 61000-4-4 | ±2kV for power supply cords ±1kV for input-output lines | ±2 kV for power supply cords Not applicable | Mains power quality should be that of a typical home or hospital environment. |
| Surge IEC 61000-4-5 | ± 0,5 kV, ± 1 kV differential mode ± 0,5 kV, ± 1 kV, ± 2 kV common mode | ± 0,5 kV, ± 1 kV differential mode Not applicable | Mains power quality should be that of a typical home or hospital environment. |
| Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11 | 0 % UT; 1 cycle at 0° 70 % UT; 25 cycles at 0° 0% UT ;250 cycles | 0 % UT; 1 cycle at 0° 70 % UT; 25 cycles at 0° 0% UT ;250 cycles | Mains power quality should be that of a typical home or hospital environment. |
| Power frequency (50/60 Hz) magnetic field IEC 61000-4-8 | 30 A/m | 30 A/m | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical hospital or home environment. |

Note: UT is the A.C. mains voltage prior to application of the test level.

Guidance and Manufacturer's Declaration - Electromagnetic Immunity: This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure that it is used in such an environment.

| IMMUNITY TEST | IEC60601 LEVEL | COMPLIANCE LEVEL | ELECTROMAGNETIC ENVIRONMENT - GUIDANCE |
|-------------------------------|---|---|--|
| Conducted RF IEC 61000-4-6 | 3V 150kHz-80MHz 6 Vrms | 3V 150kHz-80MHz 6 Vrms | Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.2 \sqrt{P}$ 150 kHz to 80 MHz $d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol:  |
| Radiated RF IEC 61000-4-3 | Amateur Radio & ISM Bands between 150kHz and 80 MHz 10V/m 80MHz-2.7GHz | Amateur Radio & ISM Bands between 150kHz and 80 MHz 10V/m 80MHz-2.7GHz | |

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from buildings, objects and people.

A: The field strength intensity of stationary transmitters such as wireless cellular/cordless phones and land mobile radios, amateur radio, AM or FM radio and television broadcasts cannot be predicted theoretically with accuracy. To evaluate the electromagnetic environment, an electromagnetic site survey should be considered. If the measured field strength in the device’s location is higher than the RF compliance level as stated in the above table, the device should be observed to verify normal operation. If abnormal device performance is observed, additional measures should be taken to relocate or reorient the device.

B: In the frequency range 150 kHz to 80MHz, the field strength should always be lower than 30 v/m.

RECOMMENDED SEPARATION DISTANCES BETWEEN PORTABLE AND MOBILE RF COMMUNICATIONS EQUIPMENT AND THIS DEVICE

This device is intended for use in electromagnetic environments in which radiated RF disturbances are controlled. The user of this device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and this device as outlined below, according to the maximum output power of the communications equipment.

| RATED MAXIMUM OUTPUT POWER OF TRANSMITTER (WATTS) | SEPARATION DISTANCE ACCORDING TO FREQUENCY OF TRANSMITTER (METERS) | | |
|---|--|---------------------------------|----------------------------------|
| | 150KHZ-80MHZ $D=1.2\sqrt{P}$ | 80MHZ-800MHZ $D=1.2\sqrt{P}$ | 800MHZ-2.5GHZ $D=2.3\sqrt{P}$ |
| 0.01 | 0.12 | 0.12 | 0.23 |
| 0.1 | 0.38 | 0.38 | 0.73 |
| 1 | 1.2 | 1.2 | 2.3 |
| 10 | 3.8 | 3.8 | 7.3 |
| 100 | 12 | 12 | 23 |


For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in Watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

APPENDIX B: CIRCUIT DIAGRAM

EN

LIMITED WARRANTY

| MODEL | WARRANTY PERIOD |
|---------|-----------------|
| OC-E80 | 3 years |
| OC-E100 | 3 years |

The manufacturer warrants that the system shall be free from defects of workmanship and materials and will perform in accordance with the product specifications for a period of 3 years from the date of sale by GCE, s.r.o. (GCE) or its affiliates to the dealer.

Filters, cannula, tubing are not covered under warranty.

Accessories are warranted to be free of defects in materials and workmanship for a period of 90 days from the time of purchase. If the product fails to perform in accordance with the product specifications, GCE, s.r.o. will repair or replace – at its discretion – the defective material or part. GCE, s.r.o. will pay customary freight charges from GCE, s.r.o. to the dealer location only. This warranty does not cover damage caused by accident, misuse, abuse, alteration, and other defects not related to material or workmanship.

GCE, s.r.o. and manufacturer disclaim all liability for economic loss, loss of profits, overhead, or consequential damages which may be claimed to arise from any sale or use of this product. Some states do not allow the exclusion or limitation of incidental or consequential damages, so the above limitation or exclusion may not apply to you.

This warranty is given in lieu of all other express or implied warranties, including the implied warranties of merchantability and fitness for a particular purpose. In addition, in no event shall the manufacturer and GCE be liable for lost profits, loss of good will, or incidental or consequential damages, even if the manufacturer and GCE have been advised of the possibility of the same. Some states or provinces do not allow the exclusion of limitation of implied warranties or the disclaimer of incidental and consequential damages. Accordingly, the laws of your state or province may give you additional protections.

To exercise your rights under this warranty, contact your local authorized GCE representative.

SysMed (CHINA) CO., Ltd.
 No.299, Baita Third Street
 Hunnan District
 110169 Shenyang
 PEOPLE'S REPUBLIC OF CHINA

EU Importer:
 GCE, s.r.o.
 Zizkova 381,
 583 01 Chotebor, Czech Republic
www.gcegroupp.com

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Shanghai International Holding Corp
GmbH (Europe), Eiffestrasse 80,
EUMEDIQ AG