



- **EXTREMELY SILENT RUN**
- **6 YEAR WARRANTY ON THE VACUUM PUMP**
- LOW MAINTENANCE REQUIREMENTS
- MAINTENANCE-FREE GERMAN VACUUM PUMP
- DESIGNED FOR CONTINUOUS RUN
- RELIABLE AND VIBRATION-FREE RUN
- VACUUM PUMP IS COOLED BY EXTRA FAN
- INNOVATIVE EUROPEAN DESIGN

- **EXTREMELY SILENT RUN**
- **6 YEAR WARRANTY ON THE VACUUM PUMP**
- LOW MAINTENANCE REQUIREMENTS
- MAINTENANCE-FREE GERMAN VACUUM PUMP
- DESIGNED FOR CONTINUOUS RUN
- RELIABLE AND VIBRATION-FREE RUN
- VACUUM PUMP IS COOLED BY EXTRA FAN
- INNOVATIVE EUROPEAN DESIGN



Overflow protection

Vacuum gauge in kPa, mmHg and Bar, colour coded

Illuminated on/off switch meets EU standard EN ISO 10079-1

Ergonomic vacuum regulator

Casing is made of very sturdy material with UV protection

Change-over valve shift the vacuum from one jar to the other with a single move – optional

Jar holder is equipped with suction circuit holder

The stainless eurorail ensures maximum flexibility. You can attach up to four jars of any size and/or any other accessories needed

MAKES LESS NOISE THAN YOUR PC!

A SPECIAL SYSTEM OF NOISE AND VIBRATION DAMPENING FOR AN EXTREMELY QUIET AND SMOOTH OPERATION WITH AN INCREDIBLY SILENT RUN

Non-marking casters with safety brakes. Antistatic casters - optional

* accessories suggestion

www.cheiron.eu

VICTORIA VERSA

The most popular high performance model of Victoria especially suitable for wards, operating rooms, recovery rooms, ICUs, obstetrics and endoscopic applications, etc.

VICTORIA LIPOS

High performance solution designed particularly for liposuction and operating rooms.



VICTORIA PORTABLE

Quiet and effective model for medical practices, wards and endoscopic applications, as well suitable for fast suction in any situation.



VICTORIA THORAX

Model specially designed for common use as well as for thoracic and abdominal surgery drainage and intensive care. The system can be used for long-term up to 30 days continuous drainage.

Thanks to a unique patented solution Victoria Thorax enables full suction output even during drainage.



OPTIONAL ACCESSORIES

Vacuum gauge in kPa, bar, cmH₂O, mmHg, colour coded



Bottle against overflow Victoria



Support ring for 2 l jar



Clamp 25x5



Easy and safe handling



Disposable secreta bag MONOKIT



0,5 – 4 l polycarbonate and polysulfone suction jars, autoclavable



MSF microbiological suction filter (20pcs/box)



Shellf



Change-over valve



IEC power outlet. It is possible to order unit with IEC power outlet



Air-operated on/off footswitch



Suction circuit holder



Storage basket



Holders allow you to place extra jars



KEYPORT Medical NZ distributor
p. 0800KEYPORT
e. sales@keyport.co.nz



CHEIRÓN a.s., Kukulova 24, Břevnov 169 00 Prague 6, Czech Republic
Branch office: Republikanska 1102/45, 312 00 Plzen, Czech Republic
Tel.: +420 377 590 433, 422, Fax: +420 377 590 435
E-mail: export@cheiron.eu, www.cheiron.eu

www.cheiron.eu

In regard of the continuous research and development Cheirón reserves the right to make changes which may affect the information contained in this brochure without prior notice.



ISO 13485