

MEDICAL HIGH PRESSURE REGULATORS

MEDIREG II®

Part No	Description
551193	Oxygen, Inlet 20,000 kpa @ 15c, Outlet x 1 SIS (400 kpa*), No flow, Swivel Yoke
551194	Nitrous Oxide, Inlet 20,000 kpa @ 15c, Outlet x 1 SIS (400 kpa*), No flow, Swivel Yoke
551195	Entonox, Inlet 20,000 kpa @ 15c, Outlet x 1 SIS (400 kpa*), No flow, Swivel Yoke
551196	Medical Air, Inlet 20,000 kpa @ 15c, Outlet x 1 SIS (400 kpa*), No flow, Swivel Yoke

*Indicates Flowing Nominal Outlet Pressure



MEDICAL HIGH PRESSURE REGULATORS

- Regulator with pressure outlet, constantly adjusted flow or with flowmeter
- Rotating pressure gauge which allows convenient reading
- Ergonomic and streamlined design
- Easy cleaning surface
- Compact and user friendly

MEDICAL HIGH PRESSURE REGULATORS

THE ADVANTAGES OF MEDIREG II®

